

8" DIAMETER TEST SIEVES

DESIGNATION		STAINLESS FRAME – STAINLESS CLOTH		BRASS FRAME – STAINLESS CLOTH			BRASS FRAME – BRASS CLOTH	
U.S.A.	METRIC	FULL HT. PART #	HALF HT. PART #	FULL HT. PART #	HALF HT. PART #		FULL HT. PART #	HALF HT. PART #
4"	00 mm	5167	—	4835	—			
3-1/2"	90 mm	5168	—	4836	—			
3"	75 mm	5169	—	4837	—			
2-1/2"	63 mm	5170	—	4838	—			
2.12"	53 mm	5171	—	4839	—			
2"	50 mm	5172	—	4840	—			
1-3/4"	45 mm	5173	—	4841	—			
1-1/2"	37.5 mm	5174	—	4842	—			
1-1/4"	31.5 mm	5175	—	4843	—			
1.06"	26.5 mm	5176	—	4844	—			
1"	25.0 mm	5177	—	4845	—			
7/8"	22.4 mm	5178	—	4846	—			
3/4"	19.0 mm	5179	—	4847	—			
5/8"	16.0 mm	5180	—	4848	—			
.530"	13.2 mm	5181	—	4849	—			
1/2"	12.5 mm	5182	—	4850	—			
7/16"	11.2 mm	5183	—	4851	—			
3/8"	9.5 mm	5184	—	4852	—			
5/16"	8.0 mm	5185	—	4853	—			
.265"	6.7 mm	5186	—	4854	—			
1/4"	6.3 mm	5187	5353	4855	5021			
No. 3-1/2	5.6 mm	5188	5354	4856	5022			
No. 4	4.75 mm	5189	5355	4857	5023			
No. 5	4.0 mm	5190	5356	4858	5024			
No. 6	3.35 mm	5191	5357	4859	5025			
*1/8"	3.15 mm	2058	3279	3281	2059			
No. 7	2.8 mm	5192	5358	4860	5026			
No. 8	2.36 mm	5193	5359	4861	5027	4529	4695	
No. 10	2.0 mm	5194	5360	4862	5028	4530	4696	
No. 12	1.7 mm	5195	5361	4863	5029	4531	4697	
No. 14	1.4 mm	5196	5362	4864	5030	4532	4698	
No. 16	1.18 mm	5197	5363	4865	5031	4533	4699	
No. 18	1.0 mm	5198	5364	4866	5032	4534	4700	
No. 20	850 micron	5199	5365	4867	5033	4535	4701	
No. 25	710 micron	5200	5366	4868	5034	4536	4702	
No. 30	600 micron	5201	5367	4869	5035	4537	4703	
No. 35	500 micron	5202	5368	4870	5036	4538	4704	
No. 40	425 micron	5203	5369	4871	5037	4539	4705	
No. 45	355 micron	5204	5370	4872	5038	4540	4706	
No. 50	300 micron	5205	5371	4873	5039	4541	4707	
No. 60	250 micron	5206	5372	4874	5040	4542	4708	
No. 70	212 micron	5207	5373	4875	5041	4543	4709	
No. 80	180 micron	5208	5374	4876	5042	4544	4710	
No. 100	150 micron	5209	5375	4877	5043	4545	4711	
No. 120	125 micron	5210	5376	4878	5044	4546	4712	
No. 140	106 micron	5211	5377	4879	5045	4547	4713	
No. 170	90 micron	5212	5378	4880	5046	4548	4714	
No. 200	75 micron	5213	5379	4881	5047	4549	4715	
No. 230	63 micron	5214	5380	4882	5048	4550	4716	
No. 270	53 micron	5215	5381	4883	5049	4551	4717	
No. 325	45 micron	5216	5382	4884	5050	4552	4718	
No. 400	38 micron	5217	5383	4885	5051			
No. 450	32 micron	5256	5389	4887	2113			
No. 500	25 micron	5218	5384	4886	5052			
No. 635	20 micron	5257	5391	2295	2296			

- Overall Height – Full: 2-5/8"; Half: 1-5/8".
- Depth to Cloth – Full: 2"; Half: 1".
- See Pages 12 and 13 for value added NIST Traceable Certification Services.

Note: 35 mesh backing cloth is available on test sieves
100 mesh and finer.

*1/8" opening is not included in ASTM E 11 Standard.

TOP COVERS		PART #
ITEM		
Stainless		8496
Stainless w. Ring		8497
Brass		8487
Brass w. Ring		8488

BOTTOM PANS		PART #
ITEM		
Stainless – Full Ht.		8491
Stainless – Half Ht.		8492
Stainless w. Skirt – Full Ht.		8493
Stainless w. Skirt – Half Ht.		8494
Brass – Full Ht.		8481
Brass – Half Ht.		8482
Brass w. Skirt – Full Ht.		8483
Brass w. Skirt – Half Ht.		8484